




Email : sales@laboratoryinstrumentindia.com

Product Name :
Users Guide ECD Kit in a Box, Spanish

Product Code :
LBNY-0023-2200016



Description :

Size: Closed A4 (Portrait) Open A3

Number of pages: Inside 68, cover 4.

Binding: Saddle stitched.

Printing: Inside 4 colors recto/verso, Cover 4 colors recto/verso plus matt lamination recto only

Paper: Inside, FSC certified 115GSM Matt Coated. Cover, FSC certified 250 GSM Matt Coated.

Technical Specification :

Users Guide ECD Kit in a Box, Spanish



Laboratory Instrument India