




Email : sales@laboratoryinstrumentindia.com

Product Name :
User Guide Recreation Kit in Box Spanish

Product Code :
LBNY-0023-220017



Description :

Paper: Inside, FSC certified 115 GSM Matt Coated. Cover, FSC certified 250 GSM Matt Coated
Printing: Inside 4 colors CMYK recto/verso, Cover 4 colors recto/verso plus matt lamination recto only
Size: Closed A4 (Portrait) Open A3
Number of pages: Inside 60, cover 4.
Binding: Saddle stitched.

Technical Specification :

Users guide Recreation Kit in a Box, Spanish



Laboratory Instrument India