




Email : sales@laboratoryinstrumentindia.com

**Product Name :**

Users guide School in a Box Kit, Spanish

**Product Code :**

LBNY-0023-2200018



**Description :**

Paper: Inside, FSC certified 115 GSM Matt Coated. Cover, FSC certified 250 GSM Matt Coated  
Printing: Inside 4 colors CMYK recto/verso, Cover 4 colors recto/verso plus matt lamination recto only  
Size: Closed A4 (Portrait) Open A3  
Number of pages: Inside 112, cover 4.  
Binding: Spiral binding (Wire O Wire).

**Technical Specification :**

Users guide School in a Box Kit, Spanish



**Laboratory Instrument India**