



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Hearing aid, adult, BTE-P, prg

**Product Code :**  
LBNY-0024-220001



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Designed to be worn behind-the-ear for people with a mild-to-severe hearing loss. Programmable behind-the-ear digital hearing aid. This electronic device amplifies and transmits sound for adults with mild-to-severe hearing loss. This device is attached to the user's ear through an earmold or dome. It can be programmed to the user's degree and type of hearing loss. Designed to provide access to speech and aid in the perception of environmental sounds.

**Technical Specification :**

Fitting range: mild-to-severe hearing loss of varying audiometric configurations  
4+ software adjustment bands  
Four additional programs  
Compression system  
Size 13 zinc-air battery  
Directional microphones  
Adaptive noise reduction  
Feedback management system  
Manual volume control  
Telecoil facility  
Battery door acts as on/off switch  
Battery door marked with the "+" symbol to indicate the location of the positive terminal  
Water-repellent coating  
Hypoallergenic casing  
Low battery alert

---

Left and right hearing aid indicators (e.g., colour coding).



**Laboratory Instrument India**