



Email : sales@laboratoryinstrumentindia.com

Product Name :
Hearing aid, adult, BTE-P, prg

Product Code :
LBNY-0024-220001



www.laboratoryinstrumentindia.com

Description :

Designed to be worn behind-the-ear for people with a mild-to-severe hearing loss. Programmable behind-the-ear digital hearing aid. This electronic device amplifies and transmits sound for adults with mild-to-severe hearing loss. This device is attached to the user's ear through an earmold or dome. It can be programmed to the user's degree and type of hearing loss. Designed to provide access to speech and aid in the perception of environmental sounds.

Technical Specification :

Fitting range: mild-to-severe hearing loss of varying audiometric configurations
4+ software adjustment bands
Four additional programs
Compression system
Size 13 zinc-air battery
Directional microphones
Adaptive noise reduction
Feedback management system
Manual volume control
Telecoil facility
Battery door acts as on/off switch
Battery door marked with the "+" symbol to indicate the location of the positive terminal
Water-repellent coating
Hypoallergenic casing
Low battery alert

Left and right hearing aid indicators (e.g., colour coding).



Laboratory Instrument India