



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Hearing aid, adult, BTE-SP, prg

**Product Code :**  
LBNY-0024-220002



www.laboratoryinstrumentindia.com

**Description :**

Digital programmable BTE (behind-the-ear) hearing aid for people with hearing loss. Designed to provide access to speech and aid in the perception of environmental sounds. This device is attached to the user's ear through an earmold. Programmable behind-the-ear digital hearing aid. This electronic device amplifies and transmits sound for adults with moderate-to-profound hearing loss. It can be programmed to the user's degree and type of hearing loss.

**INTENDED USE**

Designed to be worn behind-the-ear for people with a moderate-to-profound hearing loss.

**Technical Specification :**

Fitting range: moderate-to-profound hearing loss of varying audiometric configurations  
Four additional programs  
Compression system  
Feedback management system  
Manual volume control  
Size 675 zinc-air battery  
Directional microphones  
Adaptive noise reduction  
4+ software adjustment bands  
Water-repellent coating  
Battery door acts as on/off switch  
Battery door marked with the "+" symbol to indicate the location of the positive terminal

---

Telecoil facility  
Hypoallergenic casing  
Low battery alert  
Left and right hearing aid indicators (e.g., colour coding).



**Laboratory Instrument India**