



Email : sales@laboratoryinstrumentindia.com

Product Name :
Hearing aid, adult, BTE-SP, prg

Product Code :
LBNY-0024-220002



www.laboratoryinstrumentindia.com

Description :

Digital programmable BTE (behind-the-ear) hearing aid for people with hearing loss. Designed to provide access to speech and aid in the perception of environmental sounds. This device is attached to the user's ear through an earmold. Programmable behind-the-ear digital hearing aid. This electronic device amplifies and transmits sound for adults with moderate-to-profound hearing loss. It can be programmed to the user's degree and type of hearing loss.

INTENDED USE

Designed to be worn behind-the-ear for people with a moderate-to-profound hearing loss.

Technical Specification :

Fitting range: moderate-to-profound hearing loss of varying audiometric configurations
Four additional programs
Compression system
Feedback management system
Manual volume control
Size 675 zinc-air battery
Directional microphones
Adaptive noise reduction
4+ software adjustment bands
Water-repellent coating
Battery door acts as on/off switch
Battery door marked with the "+" symbol to indicate the location of the positive terminal

Telecoil facility
Hypoallergenic casing
Low battery alert
Left and right hearing aid indicators (e.g., colour coding).



Laboratory Instrument India