



Email : sales@laboratoryinstrumentindia.com

**Product Name :**

Hearing aid, adult, BTE-P, preprg

**Product Code :**

LBNY-0024-220003



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Designed to provide access to speech and aid in the perception of environmental sounds. This device is designed to be used in settings without appropriate resources and facilities available for hearing aid fitting and programming. Preprogrammed behind-the-ear digital hearing aid. This electronic device amplifies and transmits sound for adults with mild-to-severe hearing loss. It is designed to be an alternative on-the-spot fitting solution. This device is attached to the user's ear through an earmold or dome. There are pre-configured sound profiles available to provide an adequate match to the user's hearing loss.

**Technical Specification :**

Fixed directional microphones  
Noise reduction  
Fitting range: mild-to-severe hearing loss of common audiometric configurations  
Size 13 zinc-air battery  
Manual volume control  
Battery door acts as on/off switch  
Battery door marked with the "+" symbol to indicate the location of the positive terminal  
4+ software adjustment bands  
Compression system  
Feedback management system  
Low battery alert  
Hypoallergenic casing  
Water-repellent coating.



**Laboratory Instrument India**