



Email : sales@laboratoryinstrumentindia.com

Product Name :

Hearing aid, adult, BTE-P, preprg

Product Code :

LBNY-0024-220003



www.laboratoryinstrumentindia.com

Description :

Designed to provide access to speech and aid in the perception of environmental sounds. This device is designed to be used in settings without appropriate resources and facilities available for hearing aid fitting and programming. Preprogrammed behind-the-ear digital hearing aid. This electronic device amplifies and transmits sound for adults with mild-to-severe hearing loss. It is designed to be an alternative on-the-spot fitting solution. This device is attached to the user's ear through an earmold or dome. There are pre-configured sound profiles available to provide an adequate match to the user's hearing loss.

Technical Specification :

Fixed directional microphones
Noise reduction
Fitting range: mild-to-severe hearing loss of common audiometric configurations
Size 13 zinc-air battery
Manual volume control
Battery door acts as on/off switch
Battery door marked with the "+" symbol to indicate the location of the positive terminal
4+ software adjustment bands
Compression system
Feedback management system
Low battery alert
Hypoallergenic casing
Water-repellent coating.



Laboratory Instrument India