



Email : sales@laboratoryinstrumentindia.com

Product Name :
Hearing aid, child, BTE-P, prg

Product Code :
LBNY-0024-220004



www.laboratoryinstrumentindia.com

Description :

Digital programmable BTE (behind-the-ear) hearing aid for people with hearing loss. Designed to provide access to speech and aid in the perception of environmental sounds. This device is typically attached to the user's ear through an earmold. Programmable behind-the-ear digital hearing aid. This electronic device amplifies and transmits sound for children with mild-to-severe hearing loss. It can be programmed to the user's degree and type of hearing loss.

Technical Specification :

Directional microphones
Adaptive noise reduction
Fitting range: mild-to-severe hearing loss of varying audiometric configurations
Size 13 zinc-air battery
Feedback management system
Manual volume control
4+ software adjustment bands
Four additional programs
Compression system
Telecoil facility
Low battery alert
Battery door acts as on/off switch
Battery door marked with the "+" symbol to indicate the location of the positive terminal
Water-repellent coating
Hypoallergenic casing

Bluetooth compatibility (depending on the brand supplied)
Paediatric fitting mode available with compatible fitting software
Left and right hearing aid indicators (e.g., colour coding)
Available in three different colours.



Laboratory Instrument India