



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Hearing aid, child, BTE-P, prg

**Product Code :**  
LBNY-0024-220004



[www.laboratoryinstrumentindia.com](http://www.laboratoryinstrumentindia.com)

**Description :**

Digital programmable BTE (behind-the-ear) hearing aid for people with hearing loss. Designed to provide access to speech and aid in the perception of environmental sounds. This device is typically attached to the user's ear through an earmold. Programmable behind-the-ear digital hearing aid. This electronic device amplifies and transmits sound for children with mild-to-severe hearing loss. It can be programmed to the user's degree and type of hearing loss.

**Technical Specification :**

Directional microphones  
Adaptive noise reduction  
Fitting range: mild-to-severe hearing loss of varying audiometric configurations  
Size 13 zinc-air battery  
Feedback management system  
Manual volume control  
4+ software adjustment bands  
Four additional programs  
Compression system  
Telecoil facility  
Low battery alert  
Battery door acts as on/off switch  
Battery door marked with the "+" symbol to indicate the location of the positive terminal  
Water-repellent coating  
Hypoallergenic casing

---

Bluetooth compatibility (depending on the brand supplied)  
Paediatric fitting mode available with compatible fitting software  
Left and right hearing aid indicators (e.g., colour coding)  
Available in three different colours.



**Laboratory Instrument India**