



Email : sales@laboratoryinstrumentindia.com

Product Name :
Hearing aid, child, BTE-SP, prg

Product Code :
LBNY-0024-220005



www.laboratoryinstrumentindia.com

Description :

This device is typically attached to the user's ear through an earmold. It can be programmed to the user's degree and type of hearing loss. Programmable behind-the-ear digital hearing aid. This electronic device amplifies and transmits sound for children with moderate-to-profound hearing loss. Designed to provide access to speech and aid in the perception of environmental sounds.

Technical Specification :

Fitting range: moderate-to-profound hearing loss of varying audiometric configurations
Directional microphones
Adaptive noise reduction
Battery door marked with the "+" symbol to indicate the location of the positive terminal
Size 675 zinc-air battery
Water-repellent coating
Telecoil facility
Low battery alert
Hypoallergenic casing
4+ software adjustment bands
Manual volume control
Battery door acts as on/off switch
Left and right hearing aid indicators (e.g., colour coding)
Available in three different colours
Bluetooth compatibility (depending on the brand supplied)
Four additional programs

Compression system
Feedback management system
Paediatric fitting mode available with compatible fitting software.



Laboratory Instrument India