



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Hearing aid, child, BTE-SP, prg

**Product Code :**  
LBNY-0024-220005



www.laboratoryinstrumentindia.com

**Description :**

This device is typically attached to the user's ear through an earmold. It can be programmed to the user's degree and type of hearing loss. Programmable behind-the-ear digital hearing aid. This electronic device amplifies and transmits sound for children with moderate-to-profound hearing loss. Designed to provide access to speech and aid in the perception of environmental sounds.

**Technical Specification :**

Fitting range: moderate-to-profound hearing loss of varying audiometric configurations  
Directional microphones  
Adaptive noise reduction  
Battery door marked with the "+" symbol to indicate the location of the positive terminal  
Size 675 zinc-air battery  
Water-repellent coating  
Telecoil facility  
Low battery alert  
Hypoallergenic casing  
4+ software adjustment bands  
Manual volume control  
Battery door acts as on/off switch  
Left and right hearing aid indicators (e.g., colour coding)  
Available in three different colours  
Bluetooth compatibility (depending on the brand supplied)  
Four additional programs

---

Compression system  
Feedback management system  
Paediatric fitting mode available with compatible fitting software.



**Laboratory Instrument India**