



Email : sales@laboratoryinstrumentindia.com

Product Name :

Wheelchair, assistnt controlld,transport

Product Code :

LBNY-0024-220006



www.laboratoryinstrumentindia.com

Description :

Wheelchair, manual assistant-controlled: Transport wheelchair.

POSTURAL SUPPORT: For older CHILD and ADULT users who can sit upright and balanced without additional support.

FUNCTION AND ENVIRONMENT: Standard wheelbase and overall length. For intermittent or short-term transport in healthcare facilities.

SUPPLIED WITHOUT CUSHION.

Inappropriate prescription and fitting of wheelchairs, cushions and accessories may harm the user's health and/or function. Therefore, this product should only be purchased after consultation with the Medical Devices Unit in Supply Division.

Technical Specification :

Size and dimensions:

Backrest height: at least 430 mm

Seat width: 500 mm

Seat depth: at least 425 mm

Overall length:

Maximum 1080 mm

Frame and seating system:

Folding frame

Rear wheel position safe. In line with back post attachment to seat

Integrated tipping lever in frame for assistant to tip wheelchair to lift front castors to clear obstacles
Push handles, integrated with back post, non-slip, ergonomic grip
Seat-to-back support angle: Not adjustable, 90 degrees with bend in back post in lower thoracic area
Seat upholstery: sling
Back support upholstery: sling
Seat angle: 3 degrees

Arm supports:

Removable arm supports
Arm support design to protect clothing and cushion from rear wheels, dust, mud, etc.
Arm support pad: Approximately 70% of seat depth.
Fixed height arm support: at least 240 mm high.

Rear wheels and parking brakes:

Two small/medium rear wheels:
Diameter x width: 400 x 25 mm OR 305 x 45 mm depending on model selected
Puncture-proof rear tyres, solid polyurethane tyres
Rear wheels parking brakes, frame mounted, push-pull type.
Fixed axles

Foot supports:

Two individual flip-up and swing-away foot supports
Seat to lower leg support angle: at least 109 degrees
Height adjustable range: at least 330-410 mm

Front castors:

Two large diameter castors
Diameter x width: 200 x 50 OR 32 mm depending on model selected
Puncture-proof, solid polyurethane castor wheels

Wheelchair weight:
Total weight: 19 kg.



Laboratory Instrument India