



Email : sales@laboratoryinstrumentindia.com

Product Name :

Wheelchair, active rough terrain, adult

Product Code :

LBNY-0024-2200012



www.laboratoryinstrumentindia.com

Description :

Wheelchair, manual, active rough terrain, rigid frame:

FUNCTION AND ENVIRONMENT: For people who self-propel. Rim-propelled. For use in rough terrain in rural outdoor settings.

Relationship of seat, back-, arm- and foot support dimensions anatomically appropriate for size of wheelchair.

POSTURAL SUPPORT: For ADULT users who can sit upright and balanced without additional support, as well as those needing additional support such as adjustment of wheelchair configuration and/or addition of postural support devices. Users who do not need high back support AND/OR arm supports.

SUPPLIED WITHOUT CUSHION.

Technical Specification :

Wheelbase and overall length:

Wheelbase longer than transport wheelchair with large rear wheels.

Overall length longer than transport wheelchair with large rear wheels

Frame: Seat

Seat: solid

Seat angle: fixed

Seat height: fixed

Frame: Back support

Rigid adjustable with low profile lateral contouring and separate padded cover, adjustable forwards, backwards, angle and height using tools

Frame: Other

Rigid frame, three wheels, pushrim propelled

Push handles, non-slip, ergonomic grip

Seat unit relative to wheelbase horizontally adjustable using tools to achieve active/safe rear wheel position

Size and dimensions:

Seat width:

4 sizes: 360-480 mm in 40 mm increments

Seat depth: Adjustable

Backrest height: Adjustable

Rear wheels and parking brakes:

Rear wheels removable without tools through twist-release mechanism

Rear wheel size: Diameter x width:

660 x 35 mm

Wheel camber: 3 degrees

Pneumatic rear tyres

Rear wheels parking brakes, frame mounted, push-pull manual type

Arm supports:

Removable low profile arm rests

Arm rests protect clothing and cushion from rear wheels, dust, mud, etc. in all configurations of rear wheel position, seat angle, seat height, etc.

Foot supports:

Seat to lower leg support angle: 90 degrees

Foot supports located behind castor wheels

Foot support: two individual rigid supports

Height adjustable: 250-410 mm

Front castors

Castor size: Diameter and width:

220 x 60 mm

Solid puncture proof front castors

Wheelchair weight:

Total weight: 20.9 kg.



Laboratory Instrument India