



Email : sales@laboratoryinstrumentindia.com

**Product Name :**

Cushion,wheelchair, posture, layered

**Product Code :**

LBNY-0024-2200014



www.laboratoryinstrumentindia.com

**Description :**

Posture support cushion, contoured, with separate top layer

For CHILD or ADULT wheelchair users who require comfort and good support to maintain stability and an aligned posture; and who are at no or low risk of pressure injuries.

Cushion and cover function to reduce heat and moisture generated when the user is sitting on the cushion.

**Technical Specification :**

Characteristics:

Separate, thin top layer

Layers of cushion not glued (to be glued following fitting) OR easily separated (hook & loop laminated) to allow for modifications to optimise posture support.

Firm, non-deforming, pre-contoured layered base

Seat bone well to be custom cut by service provider on fitting OR pre-shaped seat bone well area which can be customised (depending on model selected)

Thigh area: level with no OR removable thigh supports (depending on model selected).



**Laboratory Instrument India**

