



Email : sales@laboratoryinstrumentindia.com

Product Name :

Cushion, wheelchair, posture, layered

Product Code :

LBNY-0024-2200014



www.laboratoryinstrumentindia.com

Description :

Posture support cushion, contoured, with separate top layer

For CHILD or ADULT wheelchair users who require comfort and good support to maintain stability and an aligned posture; and who are at no or low risk of pressure injuries.

Cushion and cover function to reduce heat and moisture generated when the user is sitting on the cushion.

Technical Specification :

Characteristics:

Separate, thin top layer

Layers of cushion not glued (to be glued following fitting) OR easily separated (hook & loop laminated) to allow for modifications to optimise posture support.

Firm, non-deforming, pre-contoured layered base

Seat bone well to be custom cut by service provider on fitting OR pre-shaped seat bone well area which can be customised (depending on model selected)

Thigh area: level with no OR removable thigh supports (depending on model selected).



Laboratory Instrument India

