



Email : sales@laboratoryinstrumentindia.com

Product Name :

Cushion,wheelchair,pressure, mod layered

Product Code :

LBNY-0024-2200015



www.laboratoryinstrumentindia.com

Description :

For CHILD or ADULT wheelchair users who are at risk of pressure injuries and who need good support to maintain stability and an aligned posture. Cushion and cover function to reduce heat and moisture generated when the user is sitting on the cushion.

Pressure-redistribution cushion with moderate/partial offloading or immersion/envelopment features: Contoured cushion with separate top layer.

Technical Specification :

Characteristics:

Firm, non-deforming, pre-contoured layered base

Pre-shaped seat bone well area which can be customised

Separate, thick top layer

Layers of cushion not glued (to be glued following fitting) OR easily separated (hook & loop laminated) to allow for modifications to optimise posture support.

Thigh area: level with no OR removable thigh supports (depending on model selected)



Laboratory Instrument India

