



Email : sales@laboratoryinstrumentindia.com

Product Name :
Strap pelvis/hip, wheelchair

Product Code :
LBNY-0024-2200018



www.laboratoryinstrumentindia.com

Description :

Pelvis strap:

ENVIRONMENT: For indoor and outdoor use.

FUNCTION AND POSTURE SUPPORT: For CHILD and ADULT wheelchair users to optimise posture support of the pelvis. Durable, adjustable straps with universal design to be fitted to a range of wheelchairs. Prevents the pelvis from sliding forward or rotating back (posterior tilt).

Technical Specification :

One unit with two separate sections

Length adjustable to optimise fit after fastened

Fastener: user-friendly, quick-release, buckle type system

Padded over thigh area

Skin protection through:

Design to prevent curling of edges, bunching, wrinkling, deforming or fraying and/or localised pressure points

Padding which contours to the body

Washable

Mounting length adjustable

Eyelets in strap for mounting OR mounting mechanism for wheelchair frame, depending on model selected.

Laboratory Instrument India

