



Email : sales@laboratoryinstrumentindia.com

Product Name :

Wheelchair, active rough terrain,child

Product Code :

LBNY-0024-2200021



www.laboratoryinstrumentindia.com

Description :

Wheelchair, manual, active rough terrain, rigid frame:

POSTURAL SUPPORT: For CHILD users who can sit upright and balanced without additional support, as well as those needing additional support such as adjustment of wheelchair configuration and/or addition of postural support devices. Users who do not need high back support AND/OR arm supports.

FUNCTION AND ENVIRONMENT: For people who self-propel. Rim-propelled. For use in rough terrain in rural outdoor settings.

Relationship of seat, back-, arm- and foot support dimensions anatomically appropriate for size of wheelchair.

SUPPLIED WITHOUT CUSHION.

Technical Specification :

Wheelbase and overall length:

Wheelbase longer than transport wheelchair with large rear wheels.

Overall length longer than transport wheelchair with large rear wheels

Frame: Seat

Seat:sling

Seat angle: fixed

Seat height: fixed

Frame: Back support

Tension adjustable

Frame: Other

Rigid frame, three wheels, pushrim propelled

Push handles, non-slip, ergonomic grip

Seat unit relative to wheelbase/rear wheel position not adjustable

Size and dimensions:

Seat width:

1 size: 250-300 mm adjustable through removal/addition of foam spacers on seat

Seat depth: Adjustable

Backrest height: Adjustable

Rear wheels and parking brakes:

Rear wheels removable without tools through twist-release mechanism

Rear wheel size: Diameter x width: 508 x 42 mm

Wheel camber: 3 degrees

Pneumatic rear tyres

Rear wheels parking brakes, frame mounted, push-pull manual type

Front castors

Castor size: Diameter and width: 150 x 40 mm

Solid puncture proof front castors

Clothing guards:

Rigid clothing guards

Clothing guards protect clothing and cushion from rear wheels, dust, mud, etc. in all configurations of rear wheel position, seat angle, seat height, etc.

Foot supports:

Seat to lower leg support angle: 90-95 degrees depending on seat depth

Foot supports located behind castor wheels

Foot support: one single rigid support

Height adjustable: 220-290 mm

Wheelchair weight:

Total weight: 11 kg.



Laboratory Instrument India