



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Flak Jacket, MALE, size M

**Product Code :**  
LBNY-0025-220002



#### **Description :**

##### **Soft Armour**

Equipped with detachable collar, adjustable shoulders, deployable groin protector and rear grab strap (for speedy evacuation of injured/casualties).

Fitted with front, rear, and side pockets for insertion of hard ballistic plates.

Pullover type with Velcro-fitted strap openings to enable maximum movement of the torso and multiple adjustments to ensure snug fitting.

Inner shell manufactured in suitable material for user comfort and hygiene, which also permits ease of laundering or other way to keep material in clean and hygienic state.

Outer carrier comprising an outer shell utilizing a robust fabric appropriate for the operational requirements, backcoated with polyurethane or similar material (for moisture and mildew resistance), and to include additional finishes for water repellancy. Black webbing to attach accessories.

#### **Technical Specification :**

Suitable for hot climates with high humidity.

Colour: Navy blue

Lightweight materials. Maximum weight not to exceed 3.6 kg (excluding hard insert plates, detachable collar, deployable groin protector, and any detachable accessories).

Size: Medium. Fits chest (around centre of chest) 95-104 cm. Fits waist (above belt) 86-95 cm.

---

# Laboratory Instrument India

