



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Flak Jacket,FEMALE,size S(B cup)

**Product Code :**  
LBNY-0025-220007



#### **Description :**

##### **Soft Armour**

Equipped with detachable collar, adjustable shoulders, deployable groin protector and rear extraction grab strap (for speedy evacuation of injured/ casualties), detachable comfort lining. Fitted with front, rear, and side external pockets for insertion of hard ballistic plates. Pullover type with Velcro-fitted strap openings to enable maximum movement of the torso and multiple adjustments to ensure snug fitting.

Inner shell manufactured in suitable material for user comfort and hygiene, which also permits ease of laundering or other way to keep material in clean and hygienic state, includes clear plastic ID sleeve. Outer carrier comprising an outer shell utilizing a robust fabric appropriate for the operational requirements, back coated with polyurethane or similar material (for moisture and mildew resistance), and to include additional finishes for water repellency. Black webbing to attach accessories.

#### **Technical Specification :**

Suitable for hot climates with high humidity.

Colour: Navy blue

Maximum weight not to exceed 3.6 kg (excluding hard insert plates, detachable collar, deployable groin protector, and any detachable accessories).

Size: Small (B-cup). For Bra sizes A, B, C. Fits chest (around centre of bust) 88-97cm. Fits waist (above belt) 71-80cm

S: weight 2,940 g

Weight: (excluding hard insert plates, detachable collar, deployable groin protector, and any detachable accessories).



**Laboratory Instrument India**