



Email : sales@laboratoryinstrumentindia.com

Product Name :
Helmet,ballistic resistant,level IIIA,M

Product Code :
LBNY-0025-2200016



Description :

Helmet, ballistic resistant, protection level NIJ IIIA, medium size (M).

Technical Specification :

Edge finished with a profiled neoprene strip. Including harness assembly, comfort foam pad in crown area, adjustable headband and comfortable chin strap with chin cup and single snap fastener. Shell comprising multiple layers of ballistic composite material of light weight - total weight of helmet not to exceed 1.7kg.



Laboratory Instrument India