



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Helmet,ballistic resistant,lev IIIA,L

**Product Code :**  
LBNY-0025-220017



**Description :**

Helmet, ballistic resistant, protection level NIJ IIIA, large size (L).

**Technical Specification :**

Edge finished with a profiled neoprene strip. Including harness assembly, comfort foam pad in crown area, adjustable headband and comfortable chin strap with chin cup and single snap fastener. Shell comprising multiple layers of ballistic composite material of light weight - total weight of helmet not to exceed 1.7kg.



**Laboratory Instrument India**