



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Motorcycle, agricultural type, 1-seater

**Product Code :**  
LBNY-0026-220001



**Description :**

Anti-theft device  
Clutch lever lock  
Hand protectors  
Headlight, front  
High air intake  
Brake light  
Cargo carrier, front & rear  
Rear view mirror, 2 pieces  
Side indicators, front & rear  
Speedometer  
Horn  
Lamp protection  
Lockable fuel tank lid  
Tail light  
Sump guard  
Tool kit, basic.

**Technical Specification :**

Maximum speed: 85 km/h

Engine:  
Stroke cycle: 2

---

Power: 6.2 KW at 6250 rpm  
Number of cylinders: 1  
Starter: Kick  
Fuel tank capacity: 11.00 litre (incl. 3.00 litre in reserve)  
Fuel consumption: 1.82 litre / 100 km

Brakes:  
Front: Drum  
Rear: Drum

Tyres:  
Front: 2.75-19-4PR  
Rear: 4.10-18-4PR

Transmission:  
Drive covering: Fully enclosed

Suspension:  
Front: Telescopic, 165 mm travel  
Rear: Adjustable dual spring/dampers, 130 mm travel

Colour(s): Red

Electrical (6V system):  
Ignition: Capacitor discharge (CDI)

Weight/Volume/Dimensions:  
Estimated volume: 0.86 cbm  
Estimated dimensions: 2110 x 1080 x 930 mm  
Estimated weight: 123 kg.



**Laboratory Instrument India**