



Email : sales@laboratoryinstrumentindia.com

Product Name :
Motorcycle, agricultural type, 1-seater

Product Code :
LBNY-0026-220001



Description :

Anti-theft device
Clutch lever lock
Hand protectors
Headlight, front
High air intake
Brake light
Cargo carrier, front & rear
Rear view mirror, 2 pieces
Side indicators, front & rear
Speedometer
Horn
Lamp protection
Lockable fuel tank lid
Tail light
Sump guard
Tool kit, basic.

Technical Specification :

Maximum speed: 85 km/h

Engine:
Stroke cycle: 2

Power: 6.2 KW at 6250 rpm
Number of cylinders: 1
Starter: Kick
Fuel tank capacity: 11.00 litre (incl. 3.00 litre in reserve)
Fuel consumption: 1.82 litre / 100 km

Brakes:
Front: Drum
Rear: Drum

Tyres:
Front: 2.75-19-4PR
Rear: 4.10-18-4PR

Transmission:
Drive covering: Fully enclosed

Suspension:
Front: Telescopic, 165 mm travel
Rear: Adjustable dual spring/dampers, 130 mm travel

Colour(s): Red

Electrical (6V system):
Ignition: Capacitor discharge (CDI)

Weight/Volume/Dimensions:
Estimated volume: 0.86 cbm
Estimated dimensions: 2110 x 1080 x 930 mm
Estimated weight: 123 kg.



Laboratory Instrument India