



Email : sales@laboratoryinstrumentindia.com

Product Name :
Motorcycle, on/off road type, 2-seater

Product Code :
LBNY-0026-220002



Description :

Standard equipment:
Anti-theft device
Headlight, front
High air intake
Horn
Brake light
Cargo carrier, rear
Speedometer
Tail light
Sump guard
Tool kit, basic
Rear view mirror, 2 pcs
Side indicators, front & rear.

Technical Specification :

Engine (123 cm³):
Stroke cycle: 2
Number of cylinders: 1
Power: 11 KW at 7000 rpm
Starter: Kick
Fuel consumption: 1.67 l/ 100 km
Fuel tank capacity: 9.50 l(incl. 1.00 l in reserve)

Maximum speed: 114 km/h

Brakes:

Front: Drum

Rear: Drum

Tyres:

Front: 2.75-21-4PR

Rear: 4.10-18-4PR

Transmission:

Drive covering: Open

Suspension:

Front: Telescopic, 200 mm travel

Rear: Adjustable mono spring/dampers, 155 mm travel

Electrical (6 V system):

Ignition: Capacitor discharge (CDI)

Colour(s):

White, black, blue

Weight/Volume/Dimensions:

Estimated volume: 1.03 cbm

Estimated dimensions: 2110 x 1165 x 865 mm

Estimated weight: 98 kg.



Laboratory Instrument India