



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Motorcycle, on/off road type, 2-seater

**Product Code :**  
LBNY-0026-220002



**Description :**

Standard equipment:  
Anti-theft device  
Headlight, front  
High air intake  
Horn  
Brake light  
Cargo carrier, rear  
Speedometer  
Tail light  
Sump guard  
Tool kit, basic  
Rear view mirror, 2 pcs  
Side indicators, front & rear.

**Technical Specification :**

Engine (123 cm<sup>3</sup>):  
Stroke cycle: 2  
Number of cylinders: 1  
Power: 11 KW at 7000 rpm  
Starter: Kick  
Fuel consumption: 1.67 l/ 100 km  
Fuel tank capacity: 9.50 l(incl. 1.00 l in reserve)

---

Maximum speed: 114 km/h

**Brakes:**

Front: Drum

Rear: Drum

**Tyres:**

Front: 2.75-21-4PR

Rear: 4.10-18-4PR

**Transmission:**

Drive covering: Open

**Suspension:**

Front: Telescopic, 200 mm travel

Rear: Adjustable mono spring/dampers, 155 mm travel

**Electrical (6 V system):**

Ignition: Capacitor discharge (CDI)

**Colour(s):**

White, black, blue

**Weight/Volume/Dimensions:**

Estimated volume: 1.03 cbm

Estimated dimensions: 2110 x 1165 x 865 mm

Estimated weight: 98 kg.



**Laboratory Instrument India**