



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Motorcycle, on/off road type

**Product Code :**  
LBNY-0026-220003



**Description :**

Standard equipment:

Anti-theft device  
Brake light  
Horn  
Rear view mirror, 2 pieces  
Cargo carrier, rear  
Headlight, front  
High air intake  
Sump guard  
Tail light  
Side indicators, front and rear  
Speedometer  
Tool kit, basic.

**Technical Specification :**

Engine (171 cm<sup>3</sup>):  
Stroke cycle: 2  
Number of cylinders: 1  
Power: 12.7KW at 7000 rpm  
Starter: Kick  
Fuel consumption: 1.82 litre / 100km  
Fuel tank capacity: 9.50 litre (incl. 1.00 litre in reserve)

---

Transmission:

6 gears

Drive covering: Open

Suspension:

Front: Telescopic, 200 mm travel

Rear: Adjustable mono spring/dampers, 155 mm travel

Maximum speed: 122 km/h

Electrical (6V system):

Ignition: Capacitor discharge (CDI)

Colour(s):

White, black, blue

Brakes:

Front: Drum

Rear: Drum

Tyres:

Front: 2.75-21-4PR

Rear: 4.10-18-4PR

Weight/Volume/Dimensions:

Estimated volume: 1.03 cbm

Estimated dimensions: 2110 x 1165 x 865 mm (L x H x W)

Estimated weight: 98 kg.



**Laboratory Instrument India**