



Email : sales@laboratoryinstrumentindia.com

Product Name :
Motorcycle, on/off road type

Product Code :
LBNY-0026-220003



Description :

Standard equipment:

Anti-theft device
Brake light
Horn
Rear view mirror, 2 pieces
Cargo carrier, rear
Headlight, front
High air intake
Sump guard
Tail light
Side indicators, front and rear
Speedometer
Tool kit, basic.

Technical Specification :

Engine (171 cm³):
Stroke cycle: 2
Number of cylinders: 1
Power: 12.7KW at 7000 rpm
Starter: Kick
Fuel consumption: 1.82 litre / 100km
Fuel tank capacity: 9.50 litre (incl. 1.00 litre in reserve)

Transmission:

6 gears

Drive covering: Open

Suspension:

Front: Telescopic, 200 mm travel

Rear: Adjustable mono spring/dampers, 155 mm travel

Maximum speed: 122 km/h

Electrical (6V system):

Ignition: Capacitor discharge (CDI)

Colour(s):

White, black, blue

Brakes:

Front: Drum

Rear: Drum

Tyres:

Front: 2.75-21-4PR

Rear: 4.10-18-4PR

Weight/Volume/Dimensions:

Estimated volume: 1.03 cbm

Estimated dimensions: 2110 x 1165 x 865 mm (L x H x W)

Estimated weight: 98 kg.



Laboratory Instrument India