



Email : sales@laboratoryinstrumentindia.com

Product Name :
Motorcycle, 2-seater

Product Code :
LBNY-0026-220004



Description :

Standard equipment:
Headlight, front
Brake light
Horn
Speedometer
Side indicators, front and rear
Rear view mirror, 2 pieces
Sump guard, engine protection
Tool kit, basic
Tail light.

Technical Specification :

Engine (124 cm³):
Stroke cycle: 4
Number of cylinders: 1
Power: 9.85KW at 9500 rpm
Starter: Kick
Fuel tank capacity: 7 litre (incl. 1.5 litre in reserve)
Maximum speed: 115 km/h

Suspension:
Front: Telescopic, 200 mm travel

Rear: Dual spring/dampers, 165 mm travel

Transmission:

6 gears

Drive covering: Open

Tyres:

Front: 2.75-21-4PR

Rear: 4.10-18-4PR

Brakes:

Front: Drum

Rear: Drum

Electrical (6V system):

Ignition: Capacitor discharge (CDI)

Colour(s):

Red, white

Weight/Volume/Dimensions:

Estimated volume: 0.66 cbm

Estimated dimensions: 1700 x 934 x 840 mm (L x H x W)

Estimated weight: 104 kg.



Laboratory Instrument India