



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Motorcycle, Suzuki, 2-seater

**Product Code :**  
LBNY-0026-220005



**Description :**

Standard equipment:  
Headlight, front  
Cargo carrier, rear  
Speedometer  
Tail light  
High air intake  
Horn  
Tool kit, basic  
Rear view mirror, 2 pieces  
Side indicators, front and rear.

**Technical Specification :**

Engine (123 cm<sup>3</sup>):  
Stroke cycle: 2  
Number of cylinders: 1  
Power: 7.11KW at 6000 rpm  
Starter: Kick  
Fuel tank capacity: 13 litre (incl. 2.5 litre in reserve)  
Maximum speed: 88 km/h

Transmission: 6 gears

---

Suspension:

Front: Telescopic, 165 mm travel

Rear: Adjustable spring / dampers

Colour(s):

White

Brakes:

Front: Drum

Rear: Drum

Tyres:

Front: 2.75-21-4PR

Rear: 4.10-18-4PR

Electrical (12V system):

Ignition: Electronic

Weight/Volume/Dimensions:

Estimated volume: 0.68 cbm

Estimated dimensions: 2170 x 1150 x 935 mm (L x H x W)

Estimated weight: 102 kg.



**Laboratory Instrument India**