



Email : sales@laboratoryinstrumentindia.com

Product Name :
Motorcycle, Suzuki TS185, 2-seater

Product Code :
LBNY-0026-220006



Description :

Standard equipment:
Headlight, front
Horn
Speedometer
Tail light
Tool kit, basic
Rear view mirror, 2 pieces
Side indicators, front and rear

Technical Specification :

Engine (180 cm³):
Stroke cycle: 2
Number of cylinders: 1
Power: 12.5KW at 5500 rpm
Starter: Kick
Fuel tank capacity: 7 litre (incl. 1.5 litre in reserve)
Maximum speed: 103 km/h

Tyres:
Front: 2.75-21-4PR
Rear: 4.10-18-4PR

Electrical (12V system):

Ignition: Electronic

Colour(s):

White, blue

Transmission: 5 gears

Suspension:

Front: Telescopic, 180 mm travel

Rear: Adjustable dual spring/dampers

Brakes:

Front: Drum

Rear: Drum

Weight/Volume/Dimensions:

Estimated volume: 0.68 cbm

Estimated dimensions: 2160 x 1125 x 835 mm (L x H x W)

Estimated weight: 102 kg.



Laboratory Instrument India