



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Motorcycle, Suzuki TS185, 2-seater

**Product Code :**  
LBNY-0026-220006



**Description :**

Standard equipment:  
Headlight, front  
Horn  
Speedometer  
Tail light  
Tool kit, basic  
Rear view mirror, 2 pieces  
Side indicators, front and rear

**Technical Specification :**

Engine (180 cm<sup>3</sup>):  
Stroke cycle: 2  
Number of cylinders: 1  
Power: 12.5KW at 5500 rpm  
Starter: Kick  
Fuel tank capacity: 7 litre (incl. 1.5 litre in reserve)  
Maximum speed: 103 km/h

Tyres:  
Front: 2.75-21-4PR  
Rear: 4.10-18-4PR

---

Electrical (12V system):

Ignition: Electronic

Colour(s):

White, blue

Transmission: 5 gears

Suspension:

Front: Telescopic, 180 mm travel

Rear: Adjustable dual spring/dampers

Brakes:

Front: Drum

Rear: Drum

Weight/Volume/Dimensions:

Estimated volume: 0.68 cbm

Estimated dimensions: 2160 x 1125 x 835 mm (L x H x W)

Estimated weight: 102 kg.



**Laboratory Instrument India**