



Email : sales@laboratoryinstrumentindia.com

Product Name :
Motorcycle, Honda, On/Off-Road

Product Code :
LBNY-0026-220008



Description :

Standard equipment:
Supplied with helmet

Engine (124 cm³):
Stroke cycle: 4
Number of cylinders: 2
Power: 17KW at 8500 rpm
Starter: Electric
Fuel tank capacity: 12 litres
Maximum speed: 120 km/h

Suspension:
Front: Telescopic fork with integrated shock absorbers
Rear: Mono shock suspension system

Transmission:
5 gears
Drive covering: Open

Technical Specification :

Tyres:
Front: 90/90-19M/C

Rear: 110/90-17M/C

Brakes:

Front: Single disk

Rear: Drum

Electrical (12V system):

Ignition: Electronic

Weight/Volume/Dimensions:

Estimated volume: 1.96cbm

Estimated dimensions: 2110 x 1132 x 820mm (L x H x W)

Estimated weight: 119kg.



Laboratory Instrument India