



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Motorcycle, Honda, On/Off-Road

**Product Code :**  
LBNY-0026-220008



**Description :**

Standard equipment:  
Supplied with helmet

Engine (124 cm<sup>3</sup>):  
Stroke cycle: 4  
Number of cylinders: 2  
Power: 17KW at 8500 rpm  
Starter: Electric  
Fuel tank capacity: 12 litres  
Maximum speed: 120 km/h

Suspension:  
Front: Telescopic fork with integrated shock absorbers  
Rear: Mono shock suspension system

Transmission:  
5 gears  
Drive covering: Open

**Technical Specification :**

Tyres:  
Front: 90/90-19M/C

---

Rear: 110/90-17M/C

**Brakes:**

Front: Single disk

Rear: Drum

**Electrical (12V system):**

Ignition: Electronic

**Weight/Volume/Dimensions:**

Estimated volume: 1.96cbm

Estimated dimensions: 2110 x 1132 x 820mm (L x H x W)

Estimated weight: 119kg.



**Laboratory Instrument India**