



Email : sales@laboratoryinstrumentindia.com

Product Name :

Ambulance, Hard Top, 3-door back door swing-out type, 8 seater, 4.2 L diesel

Product Code :

LBNY-0026-2200010



Description :

Standard specification for the Ambulance:

Attendant seat, for 5 persons

Heater, w/ defroster

Main stretcher

Medicine box

Power steering

Beacon lamp, blue or red, 1 pc.

Fire extinguisher

Floor, plywood covered with linoleum

Siren, motor driven

Room lamp, 10W

Sub fuel tank, 90 l

Standard equipment:

Air cleaner, cyclone, w/ pre-cleaner

Electrical system 12v

Engine under cover

Fire extinguisher, 1kg

Free wheel hub, manual

Front door fixed sash w/o vent window

Air conditioner non CFC

Assist grip, FR 2

Bumper rear, painted
Cigarette lighter
Door mirrors
Headlamp, 2 halogen
Heater, front
High altitude compensator
Lockable fuel tank lid/cap
Mud guards, front
Fuel filter for Diesel
Fuel tank 90L & extra fuel tank, 90L
Glove box w/o key
Head restraints, 2 front
Power steering
Radio, AM, w/ 1 sp.
Rear quarter glass slide
Rear seat belts NR 2x8
Rear seat parallel for 10 seats (seat belt is only for 8 seats)
Seat belts, front (ELR3-2-3)
Steering column, tilt + collapsible
Sun visor, 2pcs
Suspension, coil front, leaf rear
Seats, front, bench
Side and rear step, steel painted
Side impact beam
Spare wheel carrier back door w/lock
Spare wheel, full size
Tool kit with jack and jack handle
Tyres: 7.50R-16-8 (E) 5.50F (081D,004D)
Upholstery, vinyl
Touch-up paint
Towing eye, front & rear
Transfer 2-4 Selectors
Transfer protector (small type)
Warning triangle
Windscreen, green, laminated.

Technical Specification :

Brakes and Tyres:
Brakes Front: Disc
Brakes Rear: Drum
Tyres Front: 7.50R-16-8(E)
Tyres Rear: 7.50R-16-8(E)
Tube: Yes
Alternator: 80
Battery: 55 Ah

Engine Information:
Mount: Front Longitudal
Formation: In-line
Assembly: OHC
Displacement: 4164
Cylinders: 6
Valves per cylinder: 2

Turbo: N/A
Fuel type: Diesel
Fuel injection: N/A
Power: 96 kw at 3800 rpm
Torque: 285 Nm at 2200 rpm
Emission control: N/A
Fuel tank capacity: 180 L.



Laboratory Instrument India