



Email : sales@laboratoryinstrumentindia.com

**Product Name :**

Ambulance, Hard Top, 3-door back door swing-out type, 8 seater, 4.2 L diesel

**Product Code :**

LBNY-0026-2200010



**Description :**

Standard specification for the Ambulance:

Attendant seat, for 5 persons  
Heater, w/ defroster  
Main stretcher  
Medicine box  
Power steering  
Beacon lamp, blue or red, 1 pc.  
Fire extinguisher  
Floor, plywood covered with linoleum  
Siren, motor driven  
Room lamp, 10W  
Sub fuel tank, 90 l

Standard equipment:

Air cleaner, cyclone, w/ pre-cleaner  
Electrical system 12v  
Engine under cover  
Fire extinguisher, 1kg  
Free wheel hub, manual  
Front door fixed sash w/o vent window  
Air conditioner non CFC  
Assist grip, FR 2

---

Bumper rear, painted  
Cigarette lighter  
Door mirrors  
Headlamp, 2 halogen  
Heater, front  
High altitude compensator  
Lockable fuel tank lid/cap  
Mud guards, front  
Fuel filter for Diesel  
Fuel tank 90L & extra fuel tank, 90L  
Glove box w/o key  
Head restraints, 2 front  
Power steering  
Radio, AM, w/ 1 sp.  
Rear quarter glass slide  
Rear seat belts NR 2x8  
Rear seat parallel for 10 seats (seat belt is only for 8 seats)  
Seat belts, front (ELR3-2-3)  
Steering column, tilt + collapsible  
Sun visor, 2pcs  
Suspension, coil front, leaf rear  
Seats, front, bench  
Side and rear step, steel painted  
Side impact beam  
Spare wheel carrier back door w/lock  
Spare wheel, full size  
Tool kit with jack and jack handle  
Tyres: 7.50R-16-8 (E) 5.50F (081D,004D)  
Upholstery, vinyl  
Touch-up paint  
Towing eye, front & rear  
Transfer 2-4 Selectors  
Transfer protector (small type)  
Warning triangle  
Windscreen, green, laminated.

---

**Technical Specification :**

Brakes and Tyres:  
Brakes Front: Disc  
Brakes Rear: Drum  
Tyres Front: 7.50R-16-8(E)  
Tyres Rear: 7.50R-16-8(E)  
Tube: Yes  
Alternator: 80  
Battery: 55 Ah

Engine Information:  
Mount: Front Longitudal  
Formation: In-line  
Assembly: OHC  
Displacement: 4164  
Cylinders: 6  
Valves per cylinder: 2

---

Turbo: N/A  
Fuel type: Diesel  
Fuel injection: N/A  
Power: 96 kw at 3800 rpm  
Torque: 285 Nm at 2200 rpm  
Emission control: N/A  
Fuel tank capacity: 180 L.



**Laboratory Instrument India**