



Email : sales@laboratoryinstrumentindia.com

Product Name :

Ambulance, 4 door w/1 sliding side door, back door lift-up type, 7-seater, 2.4 L diesel engine

Product Code :

LBNY-0026-220011



Description :

Standard specification for the Ambulance:

Main stretcher with caster

Room lamp, 1 pc

Attendant seat, 5 persons

Electric Siren Amplifier & Speaker, 1 tone

Fire extinguisher

Beacon lamp, red

Floor, plywood & linoleum covered floor.

Standard Equipment:

Air cleaner, cyclone type

Clock, built-in trip meter

Electrical system, 12V

Fire extinguisher, 1kg install

Front headrests, 3pcs for LHD, 2pcs for RHD

Air conditioner, front and rear (Non CFC)

Assist grips, front:3

Battery, 85D26Rx2(55AHx2)

Beacon lamp(red, front 55W)

Cigarette lighter

Head lamp, sculptured 2

Heater, front & rear

Inside rear view mirror, day and night
Intermittent wiper, intermittent
Front seat: 3persons
Front seat belt, 3P ELRx2 + NRx1
Front wheel brake, disc (15 inches)
Fuel lid opener, cable type
Plywood & linoleum covered floor
Power steering
Pre-air cleaner
Radio, AM/FM, w/ 2 sp.
Lockable fuel tank
Main stretcher w/ caster
Mud guards, w/o
Outside rear view mirror, 2pcs
Rear door inside handle
Rear seat, w/o
Rear wheel brake, drum
Room lamp(8W) : front + center + rear
Seat material, fabric
Steering wheel, 2-spoke urethane
Sub stretcher, folding
Sun visors, 2 pcs.
Suspension, fr:double wishborn, rr:rigid leaf
Tool kit with jack & Jack handle
Touch-up paint
Seats, attendant, folding type, 4 persons
Side door beam, front and slide door
Siren amplifier, electric & speaker (1 tone type)
Spare wheel, 1pc (full size)
Steering column, manual tilt
Towing eye, front and rear
Wheel cap, full cap
Wheel house cover
Towing hook, fr: 1
Transport hook, fr: 2
Triangle caution plate, 1pc loaded
Tyres: 195R15C 8PR Tubeless 15x6JJ
Windshield glass, green (adhesion).

Technical Specification :

Brakes and Tyres:
Brakes Front: Disc
Brakes Rear: Drum
Tube: No
Alternator: 120
Battery: 60 Ah
Tyres Front: 195R15C 8PR
Tyres Rear: 195R15C 8PR

Engine Information:
Mount: Front Longitudal
Displacement: 2986
Assembly: OHC

Power: 68 kw at 4000 rpm
Torque: 197 Nm at 2400 rpm
Cylinders: 4
Valves per cylinder: 2
Formation: In-line
Turbo: N/A
Emission control: N/A
Fuel tank capacity: 70 L
Fuel type: Diesel
Fuel injection: EFI.



Laboratory Instrument India