




Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Automatic Gas cutting M/C

**Product Code :**  
LBNY-0004-820061



**Description :**

Automatic Gas cutting M/C

**Technical Specification :**

Automatic Gas cutting M/C

The main body is made of aluminum by die-casting minutely.

Speed system is adjusted by integrated circuit, keeping running steadily.

Oxygen acetylene hose is made of high quality rubber for durable use.

Various kinds of spare parts can be cut as per pattern.

High precision in cutting with well smooth surface.

Especially for batch producing of components in same pattern.

Applied to shipbuilding, oil production, metallurgy, metal structure and other manufacturers etc.

**Specifications:**

Lightning input voltage: 220V,50-60Hz

Cut thickness: 6-50 mm

Improvement angle: 0-45°

Cutting speed:100-1000 mm/min

Round cut diameter: 360-1200 mm

**Accessory**

Cutting crater

Cutting workbench

---

Welding Mask



**Laboratory Instrument India**