



Email : sales@laboratoryinstrumentindia.com

Product Name :
2WD,diesel,15 seater

Product Code :
LBNY-0026-2200012



Description :

Standard equipment:
Air cleaner, cyclone type
Air conditioner, front and rear (Non CFC)
Assist grips, front:3
Battery, 85D26Rx2(55AHx2)
Front headrests, 3pcs for LHD, 2pcs for RHD
Front seat: 3persons
Front seat belt, 3P ELRx2 + NRx1
Front wheel brake, disc (15 inches)
Beacon lamp (red, front 55W)
Cigarette lighter
Clock, built-in trip meter
Electrical system, 12V
Fire extinguisher, 1kg install
Lockable fuel tank
Main stretcher w/ caster
Mud guards, w/o
Outside rear view mirror, 2pcs
Fuel lid opener, cable type
Head lamp, sculptured 2
Heater, front & rear
Inside rear view mirror, day and night
Intermittent wiper, intermittent

Rear seat, w/o
Rear wheel brake, drum
Room lamp (8W): front + center + rear
Seat material, fabric
Seats, attendant, folding type, 4 persons
Plywood & linoleum covered floor
Power steering
Pre-air cleaner
Radio, AM/FM, w/ 2 sp.
Rear door inside handle
Sub stretcher, folding
Sun visors, 2 pcs.
Suspension, fr:double wishborn, rr:rigid leaf
Tool kit with jack & Jack handle
Touch-up paint
Side door beam, front and slide door
Siren amplifier, electric & speaker (1 tone type)
Spare wheel, 1pc (full size)
Steering column, manual tilt
Steering wheel, 2-spoke urethane
Tyres: 195R15C 8PR Tubeless 15x6JJ
Wheel cap, full cap
Wheel house cover
Towing eye, front & rear
Towing hook, fr: 1
Transport hook, fr: 2
Triangle caution plate, 1pc loaded
Windshield glass, green (adhesion).

Technical Specification :

Engine Information:
Mount: Front Longitudal
Valves per cylinder: 2
Formation: In-line
Displacement: 2986
Cylinders: 4
Assembly: OHC
Turbo: N/A
Fuel type: Diesel
Fuel injection: EFI
Power: 68 kw at 4000 rpm
Torque: 197 Nm at 2400 rpm
Emission control: N/A
Fuel tank capacity: 70 L

Brakes and Tyres:
Brakes Front: Disc
Brakes Rear: Drum
Alternator: 120
Battery: 60 Ah
Tyres Front: 195R15C 8PR
Tyres Rear: 195R15C 8PR
Tube: No



Laboratory Instrument India