



Email : sales@laboratoryinstrumentindia.com

Product Name :

Bicycle,heavy duty,roadster,gents

Product Code :

LBNY-0026-220014



Description :

Frame size: 560mm/22".
Wheel size: 710mm/28".
Tyre driven dynamo lighting set
Partly covered chain case
Heavy duty front and rear carriers
Imitation leather saddle
Double stand
Reflectors on pedals
Tyre inflator
Coaster brakes
Bell
Frame lock
Tool bag with tools
Colours: Blue or black.

Technical Specification :

Weight/Volume/Dimensions:
Estimated volume: 0.190cbm
Estimated weight: 26.8kg



Laboratory Instrument India