




Email : sales@laboratoryinstrumentindia.com

Product Name :
Truck,pallet,manual,2000kg

Product Code :
LBNY-0027-220001



Description :

Manually-operated pallet truck, with hydraulically operated lever action and a lifting capacity of 2000 kg. The truck has a 3-point control trigger; bogie mounted vulcanized wheels and a strong vinyl-coated towing/steering handle.

Maximum lifting height: 205 mm.

Maximum fork width: 520 mm.

Length of forks: 1125 mm.

Technical Specification :

Weight/Volume/Dimensions:

Estimated volume: 0.03 cbm

Estimated dimensions: 1.13 x 0.53 x 0.50 m

Estimated weight: 77 kg



Laboratory Instrument India