



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Barrow, trolley, sack

**Product Code :**  
LBNY-0027-220002



**Description :**

The frame is painted, made of tubular steel construction with plastic grip handles. A two wheeled sack-barrow of maximum capacity 200 kg.

The barrow is fitted with solid rubber wheels running on ball bearings, size 200 x 50 mm.

Barrow dimensions: width, 520 mm; toe plate, 350 x 200 mm; height, 1,070 mm.

**Technical Specification :**

Weight/Volume/Dimensions:

Estimated volume: 0.003 cbm

Estimated dimensions: 1,100 x 600 x 350 mm

Estimated weight: 22 kg.



**Laboratory Instrument India**