



Email : sales@laboratoryinstrumentindia.com

Product Name :
Sling,chain,3000kg capacity

Product Code :
LBNY-0027-220003



Description :

Max. lifting capacity 3000 kg, length of each chain leg 1000 mm.
A double leg type chain sling, 10 mm diameter with alloy steel links, complete with pear link and two slip hooks.

Technical Specification :

Weight/Volume/Dimensions:
Estimated volume: 0.002 cbm
Estimated dimensions: 250 x 100 x 100 mm
Estimated weight: 16.5 kg.



Laboratory Instrument India