



Email : sales@laboratoryinstrumentindia.com

**Product Name :**

Warehouse pref.,portable,GI/Alu+PVC,24x10m

**Product Code :**

LBNY-0027-220004



**Description :**

The Warehouse has one large hooded ventilation panel at each gable end fitted with mosquito netting. A prefabricated, (prefab) portable kit made from galvanized steel or aluminium frame and covered with strong white PVC coated Polyester sheeting.

**Technical Specification :**

The concept of the warehouse is modular (usually with 4.0m modules) of 24-meter length and 10 meter width, designed for easy and rapid erection. The modular concept allows for extending the length if needed. The height of the vertical side walls is a minimum 3.0m and the height at the apex is about 4.0m.

**The frame:**

The frame is made of aluminium or steel sections. Steel is hot-dipped galvanized after forming and welding. Galvanization is in accordance with approved international standards. The frame should be designed to withstand winds of up to 35m/second and a snow load of minimum 50 kp/m<sup>2</sup>.

**Cover sheeting:**

PVC coated polyester fabric sheeting, UV stabilized, self extinguishing, waterproof, suitable for hot and cold climates resisting temperatures of -30 to +70 °C.

**Doors:**

The Warehouse is fitted with two doors; one fold-up door and one fold-away door made of PVC, minimum size 4 x 4meters.



**Laboratory Instrument India**