



Email : [sales@laboratoryinstrumentindia.com](mailto:sales@laboratoryinstrumentindia.com)

**Product Name :**  
Colorimeter dgt pocket,0-2mg/l Chlorine

**Product Code :**  
LBNY-0028-220002



**Description :**

Light weight Pocket Colorimeter.

It provides Liquid Crystal Digital Display and direct concentration read-out.

Battery operated using four (4) AAA alkaline batteries.

A hand-held filter photometer factory programmed with the calibration for a single parameter (Chlorine), using the DPD method.

Range of chlorine analysis: 0-2mg/litre chlorine.

(Free and total).

Reagents are needed for the tests.

Smallest increment 0.01mg/litre.

**Technical Specification :**

Weight/Volume/Dimensions:

Estimated volume: 270cdm

Estimated weight: 0.2kg.



**Laboratory Instrument India**

