



Email : sales@laboratoryinstrumentindia.com

Product Name :
Spot Welder

Product Code :
LBNY-0004-820062



Description :

Spot Welder

Technical Specification :

Spot Welder

Specifications:

Rated capacity(KVA): 50

Rated primary voltage (V): 220/380/44

Max welding input (KVA): 150(129)

Short-circuit current(A):20,000

Welding current(A): 18,000-20000

Rated pressing force(kgf): 600-700

Throat depth(mm): 500-550

Neck spacing(mm):210

Allowable usage rate(%): 5.5

Electrode stroke(Hz): 50

Electrode mirror(V): 16

Electrode hold diameter(\varnothing): 25-30

Electrode hold length(mm): i80-200

Horn diameter(\varnothing): 60-70 mm

Accessory

Electrode holder

Electrode tip
Welding Mask



Laboratory Instrument India