



Email : sales@laboratoryinstrumentindia.com

Product Name :
Bags,Water Sampling,Sterile,Stand,100ml

Product Code :
LBNY-0028-220003



Description :

Water sampling bag, polyethylene, sterile, self-stand-up bag type. The bag has a leak-proof closure with folding tap closing, 100ml and 4oz. Fill lines and a write-on strip.

Water Sampling bags, 100ml capacity, Sterile, Stand-up, easy and safe closing, without Thiosulphate tablets.

Technical Specification :

Dimensions

3" Wx7-1/4"L (7.5x18.5cm)

2.5mil(0.064mm)thick,

4oz(118ml)capacity.



Laboratory Instrument India