



Email : sales@laboratoryinstrumentindia.com

**Product Name :**  
Bact.testing,E-Coli,Compact Dry Plates

**Product Code :**  
LBNY-0028-220004



**Description :**

Compact Dry E. coli/Coliform Count (EC) is a Direct Culture based ready-to-use test method for Rapid Enumeration of Escherichia coli and coliform bacteria.

Compact Dry E. coli/Coliform Count (EC) plates for bacteriological testing of drinking water.

The medium should be rehydrated with 1mL diluted sample inoculated onto the center of the self-diffusible medium, allowing the solution to diffuse by capillary action. The plate can be incubated at 35[degrees]C for 20-24h and the colonies counted without any further working steps.

The plates are pre-sterilized, contain culture medium and a cold water-soluble gelling agent

Medium contains 2 kinds of chromogenic enzyme substrates that develop red/pink color for Coliform and white/blue color for E.coli.

Storage: Keep at room temperature (1-30 degrees C). Refrigerated storage is not required.

**Technical Specification :**

Shelf life: 18 months after manufacturing.



**Laboratory Instrument India**

